

ODCZYNNIKI ANTYGLOBULINOWE

Produkt	Nr kat.	Opakowanie	Charakterystyka
Anti-IgG,-C3d	5441013 5441023	10 ml 10 x 10 ml	<i>NovaClone</i> ™ Anti-IgG :5H4+8D2-8 Anti-C3c: 86 5A2 Anti-C3d :139 4B4
Anti-IgG,-C3d (green)	5451013 5451023	10 ml 10 x 10 ml	j.w., zawiera zielony barwnik (acid blue/acid yellow)
Anti-IgG	5461013 5461023	10 ml 10 x 10 ml	<i>NovaClone</i> ™ Anti-IgG : 5H4+8D2-8
Anti-IgG (green)	5470013 5470023	10 ml 10 x 10 ml	j.w., zawiera zielony barwnik (acid blue/acid yellow)
Anti-IgA	0003002	2 ml	odczynnik poliklonalny
Anti-IgM	0003003	2 ml	odczynnik poliklonalny
Anti-C3d	5491053	3 ml	<i>NovaClone</i> ™ Anti-C3d :139 4B4